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Type of Action: Correction

Month DD, YYYY (to be adapted locally)

Dear Dialysis Provider: (to be adapted locally)

This communication is a follow-up to the attached Correction notification you received from Vantive, formerly Baxter, dated October 16, 2024 (to be adapted). Vantive issued the previous communication to inform patients about potential safety risks related to exposure to non-dioxin-like (NDL) polychlorinated biphenyl acids (PCBAs) and NDL polychlorinated biphenyls (PCBs) when using certain peritoneal and hemodialysis devices. This initial notice was based on publicly available information at that time. At that time Vantive was in the process of evaluating whether these same risks are present with the AK 98 hemodialysis machines (“AK 98”).

Vantive now has a Toxicological Risk Assessment (TRA) with respect to the actual relative risks associated with the use of any AK 98 machines containing peroxide-cured silicone tubing components. These test results currently demonstrate that the presence of the peroxide-cured silicone tubing components is unlikely to pose a risk for patients with a body weight greater than 25kg, irrespective of the run-time of the machine (note: 25kg is the minimum patient weight indicated for AK 98). Accordingly, Vantive has determined that no corrections are needed for all devices currently with customers, and patients may continue to perform therapy with the current configuration of AK 98.

Affected Product (to be adapted)

| Product Code | Product Description | Serial Numbers |
|--------------|----------------------------|---------------------------------|
| 115244 | AK 98, 230V, EFFICIENT | All serials lower than AK078305 |
| 115248 | AK 98, 230V, BIO | All serials lower than AK078305 |
| 115249 | AK 98, 115V, BIO | All serials lower than AK078305 |
| 115250 | AK 98, 230V, SELF-CARE | All serials lower than AK078305 |
| 955403 | AK 98 V2 230V | All serials lower than AK078305 |
| 955404 | AK 98 V2 230V, SELF-CARE | All serials lower than AK078305 |
| 955406 | AK 98 V2 115V, BIO VERSION | All serials lower than AK078305 |

| | | |
|--------|--------------------------|---------------------------------|
| 955603 | AK 98V3 230V BIO VERSION | All serials lower than AK078305 |
| 955604 | AK 98V3 230V SELF-CARE | All serials lower than AK078305 |
| 955989 | AK 98 AK98 V3RU | All serials lower than AK078305 |

Hazard Involved

Vantive’s latest test results demonstrate that the presence of the peroxide-cured silicone tubing components in AK 98 is unlikely to pose a risk for patients with a body weight greater than 25kg, irrespective of the run-time of the machine. To date, Vantive has not received any complaints determined to be related to this issue.

Actions to be Taken by Customers

1. No further action is necessary; all AK 98 machines, regardless of run-time, may continue to be used and no corrections are required.
2. Complete the enclosed customer reply form and return it to Vantive by either scanning and e-mailing it to (insert local contact information) or sending it by post to (insert local contact information), even if you don’t have any inventory. Returning the customer reply form promptly will confirm your receipt of this notification and prevent you from receiving repeat notices.
1. If you purchased this product from a distributor, please note that responding with the Vantive customer reply form is not applicable. If a response is requested by your distributor or wholesaler, please respond to them according to their instructions.
2. If you distributed this product to other facilities or departments within your institution, please forward a copy of this communication to them.
3. If you are a dealer, wholesaler, distributor/reseller, or original equipment manufacturer (OEM) that distributed any affected product to other facilities, please notify your customers of this notification in accordance with your customary procedures. (to be adapted locally)

Further Information and Support

For general questions regarding this communication or any product issue you are experiencing, contact Vantive at (insert local contact information), between the hours of (insert local information).

The local Ministry of Health (MOH) has been notified of this action. (to be adapted locally)

We apologize for any inconvenience this may cause you and your staff.

Sincerely,

Amy McKernan (to be adapted)
Senior Director, Product Quality (to be adapted)
Vantive US Healthcare LLC (to be adapted)

Enclosures: Vantive Customer Reply Form
Dialysis Provider letter dated October 16, 2024 (to be adapted)
Home Patient Letter (to be adapted)

(Customer communication)

CUSTOMER REPLY FORM related to Product Correction letter dated XXXXXX (to be completed locally)
Product Name: AK 98 (to be adapted locally)
Product code: (to be adapted locally)
Serial Number: All

Please complete and return one copy of this form per facility by e-mail (_____) as confirmation that you have received this notification. *(Can be adapted locally).*

| | |
|---|--|
| Facility Name and Address: | |
| Reply Confirmation Completed By (Please Print): | |
| Title (Please print): | |
| Email and/or Telephone Number (including Area Code): | |

Your signature below indicates that you have received the attached letter; performed the actions as outlined in the letter as needed; and disseminated this information to staff and other services or facilities as applicable.

| | |
|--|-------|
| Signature/Date: REQUIRED FIELD | <hr/> |
|--|-------|