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MiniCap Extended Life PD transfer

FA Number: FAV-2024-007

Manufacturer: Baxter Healthcare SA (CH-MF-000026124)

Type of Action: Correction

Month DD, YYYY (to be adapted locally)

Dear Peritoneal Dialysis Patient:

Baxter Healthcare Corporation (Baxter) is issuing a Correction for the MiniCap Extended Life PD transfer sets listed below, which are manufactured with peroxide-cured silicone tubing that transfers solution to and from the patient.


Baxter is aware of several recalls by other manufacturers related to the potential risk of exposure to non-dioxin-like (NDL) polychlorinated biphenyl acids (PCBAs) and NDL polychlorinated biphenyls (PCBs) when using certain peritoneal dialysis and hemodialysis devices. The source of the NDL PCBAs and/or NDL PCBs in those recalls was due to the manufacturing process of the silicone tubing.

Baxter is in the process of evaluating whether these same risks are present with its PD transfer sets. At this time Baxter does not have data to definitively conclude whether there is a safety risk. Therefore, Baxter is informing you of the potential patient safety risk while our evaluation is in process.

While this evaluation is ongoing, Baxter is also in the process of transitioning certain product codes of the PD transfer sets from peroxide-cured silicone tubing to platinum-cured silicone tubing. Available information indicates that NDL PCBAs and NDL PCBs are not detected in medical devices with this modified version of silicone tubing. Please note that the exact timing of this transition will vary by geographic region, and that Baxter will continue to make the existing peroxide-cured silicone tubing configuration of the transfer sets available in your country until this transition occurs, as there is currently no definitive data to demonstrate that a patient safety risk exists. The purpose of this letter is to inform you of this current status, and to let you know that as additional data becomes available, Baxter will provide you with further communication if any mitigation measures are necessary for transfer sets with the peroxide-cured silicone tubing components.

Affected Product (to be adapted)

Product Code	Product Description	Lot Numbers
5C4482	Transfer Set (MiniCap Extended Life PD Transfer Set with Twist Clamp)	All lots within expiry
R5C4482	Transfer Set (MiniCap Extended Life PD Transfer Set with Twist Clamp)	
R5C4482E	Transfer Set (MiniCap Extended Life PD Transfer Set with Twist Clamp)	
R5C4483	Transfer Set (MiniCap Extended Life PD Transfer Set with Twist Clamp)	

R5C4484	Transfer Set (MiniCap Extended Life PD Transfer Set with Twist Clamp)	
---------	---	--

Hazard Involved

Polychlorinated biphenyls are persistent organic pollutants that have a negative impact on the ecosystem and all living beings and continue to represent a serious risk to human health. The risks include neuropsychological, neurobehavioral deficits, dementia, immune system dysfunctions, cardiovascular diseases, cancer and harmful effects on the reproductive system. Baxter is in the process of evaluating whether these risks are present with the MiniCap Extended Life PD transfer sets. At this time Baxter does not have data to definitively conclude whether there is a safety risk. To date, Baxter has not received any complaints related to this issue.

Actions to be Taken by Patients

1. **Please continue your prescribed PD therapy.** If you have any questions about your PD therapy, please contact your doctor and/or nurse.
2. Complete the enclosed customer reply form and return it to Baxter by either scanning and e-mailing it to (insert local contact information) or sending it by post to (insert local contact information), even if you don't have any inventory. Returning the customer reply form promptly will confirm your receipt of this notification and prevent you from receiving repeat notices.

Further Information and Support

If you have any questions about your PD therapy, please contact your doctor and/or nurse. (to be adapted)

For general questions regarding this communication or any product issue you are experiencing, contact Baxter at (insert local contact information), between the hours of (insert local information).

We apologize for any inconvenience this may cause you and your staff.

Sincerely,

Name (to be adapted locally)

Title (to be adapted locally)

Baxter Healthcare Corporation (to be adapted locally)

Enclosures: Baxter Customer Reply Form

Confirmation of receipt of communication

(CORRECTION LETTER DATED XX (TO BE COMPLETED LOCALLY))

DEVICE NAME: EXTENDED LIFE PD TRANSFER SET, MINICAP EXTEND LIFE

Product code : 5C4482, R5C4482, R5C4482E, R5C4483, R5C4484 (to be completed locally)

Batch numbers : All lots within expiry

Please complete and return one copy of this form by e-mail (_____) as confirmation that you have received this notification.
(Can be adapted locally)

Name:	
Email and/or Telephone Number (including Area Code):	

Signature/Date: REQUIRED FIELD	_____
-----------------------------------	-------

I have received the above-mentioned letter and performed the actions outlined in the letter.